

Recall Information

Medline Bed Assist Bar

Recall Date: 5/30/2024

SKU #	Internet #	Product Description	UPC	Model #
1001096203	204993897	Medline Bed Assist Bar	884389102865	BA6800



Name of Product: Medline Bed Assist Bar Adult Bed Rails

Hazard: When the recalled bed rails are attached to an adult's bed, users can become entrapped within the bed rail or between the bed rail and the side of the mattress. This poses a serious entrapment hazard and risk of death by asphyxiation.

Remedy: Consumers should immediately stop using the recalled bed rails and contact Medline for a refund.

Consumer Contact: Medline toll-free at 866-359-1704 from 8 a.m. to 5 p.m. CT Monday through Friday, by email at recalls@medline.com, or online at <https://athome.medline.com> and click on "Product Recalls" at the bottom of the page for more information.